

## Quantification of peptide adsorption on lipid nanocapsules without prior purification or separation process

Claire Gazaille, Marthe Akiki, Joël Eyer, Guillaume Bastiat

### ► To cite this version:

Claire Gazaille, Marthe Akiki, Joël Eyer, Guillaume Bastiat. Quantification of peptide adsorption on lipid nanocapsules without prior purification or separation process. Journée scientifique de la SFR ICAT 4208, Jun 2019, Angers, France. , 2019. hal-02616080

### HAL Id: hal-02616080 https://univ-angers.hal.science/hal-02616080

Submitted on 24 May 2020

**HAL** is a multi-disciplinary open access archive for the deposit and dissemination of scientific research documents, whether they are published or not. The documents may come from teaching and research institutions in France or abroad, or from public or private research centers. L'archive ouverte pluridisciplinaire **HAL**, est destinée au dépôt et à la diffusion de documents scientifiques de niveau recherche, publiés ou non, émanant des établissements d'enseignement et de recherche français ou étrangers, des laboratoires publics ou privés.





Institut national de la santé et de la recherche médicale





# Quantification of peptide adsorption on lipid nanocapsules

## without prior purification or separation process

Gazaille C, Akiki M, Eyer J, Bastiat G\* Translational Micro and Nanomedicines (MINT, INSERM 1066 / CNRS 6021), UNIV ANGERS, Angers, France

\*corresponding author: guillaume.bastiat@univ-angers.fr

To quantify the encapsulation rates of active agents (drug, protein, etc.) in nanoparticles, the separation of non-loaded agents from the loaded ones is often reported in literature, using dialysis or centrifugal filters.<sup>1,2</sup> However, these techniques are more and more questioned, and the low quality of the separation could lead to over or underestimated data.<sup>3,4</sup> To avoid the biases, a method was developed to quantify free peptide proportions without prior separation step from peptide adsorbed at the nanoparticle surface, using NFL-TBS.40-63 (NFL) and lipid nanocapsules (LNCs) as peptide and nanoparticle models, respectively.



(3.5 minutes) and NFL (7.7 minutes) on size exclusion column.





Thanks to a size exclusion method, the separation between nanoparticles and active agent and their quantification are performed simultaneously. This single step is a real benefit for analysis accuracy and robustness, minimizing the biases compared to sequential steps for separation and quantification. In this particular case, the best combinations of concentrations/types of LNCs and NFL were determined to achieve a total NFL adsorption at LNC surface. This technique could be a promising tool to characterize other nanoparticle/active agent interactions.

## References

1. Chen, L. et al. Int. J. Biol. Macromol. 2019, 254–261. 2. Abriata, J. P. et al. Mater. Sci. Eng. C.2019, 347–355. 3. Johnsen, E. et al. J. Pharm. Biomed. Anal. 2016, 106–111. 4. Garvin, A. et al. J. Forensic Sci. 2013, S173–S175.

5. Heurtault, B. et al. Pharm Res. 2002, 875-80. 6. Heurtault, B. et al. Biomaterials. 2003, 4283–4300. 7. Berges, R. et al. PLoS One. 2012, e49436.

## Acknowledgments

Reconnue d'utilité publique

